

Cut Sets for TOP L

71850 - TC-1 Sergnano

Cut set	Order	Cut set event descriptions	Frequency (occasioni/anno)	Unavailability
EV121.EV125	2	Mancato avvio ventilatori E-101 TSHH151B intervento logica	7,943e-4	1,085e-5
EV121.EV166	2	Mancato avvio ventilatori E-101 Mancata apertura HSV107	6,649e-4	1,67e-4
EV123.EV166	2	Logica Mancata apertura HSV107	4,694e-4	6,412e-6
EV121.EV124	2	Mancato avvio ventilatori E-101 TE151B	4,663e-4	6,37e-6
EV121.EV165	2	Mancato avvio ventilatori E-101 Mancata chiusura HSV106	3,172e-4	7,959e-5
EV122.EV166	2	TE151 Mancata apertura HSV107	2,756e-4	3,765e-6
EV123.EV165	2	Logica Mancata chiusura HSV106	2,238e-4	3,056e-6
EV122.EV165	2	TE151 Mancata chiusura HSV106	1,314e-4	1,794e-6
EV123.EV125	2	Logica TSHH151B intervento logica	5,935e-5	4,166e-7
EV122.EV125	2	TE151 TSHH151B intervento logica	3,484e-5	2,446e-7
EV123.EV124	2	Logica TE151B	3,484e-5	2,446e-7
EV122.EV124	2	TE151 TE151B	2,046e-5	1,436e-7